



# **CORRECTIONAL DENTAL ASSOCIATES**

## **Correctional Dental Prosthetic Tracking & Compliance System Installation Manual**

**The Correctional Dental Prosthetic Tracking & Compliance System (CDPTCS) requires a minimum of a Pentium IV, 256 megabytes RAM, 30 megabytes free hard disk space, and a Windows 2000 or Windows XP operating system.**

### **Step 1**

Install CDPTCS Installation Disk in CD Rom drive of computer.

### **Step 2**

Click Start, then Click Run

In the Run Dialog Box type "D: setup.exe", Click "OK"

The CDPTCS Install Window will appear and the install program will start copying the installation files to the computer.

CDPTCS Setup Window appears

Click “OK” to Continue. Next click installation button on left (button icon looks like a computer with a software box) to install program in default folder C:\Program Files\CDPTCS\

(To change the directory click “Change Directory” button.)

### **Step 3**

Click Continue at the CDPTCS – Choose Program Group Window

The CDPTCS Setup and database installation will begin.

*(Click “Ignore” on any error messages that may come up in the section of the installation.)*

Click ‘OK” Button to “CDPTCS Setup was completed successfully”

### **Step 4**

Click **Start** then click **All Programs**

Go to the CDPTCS folder and click on the CDPTCS application Icon

Click Start > Control Panel > Administrative Tools

Double Click **Data Sources ODBC**

ODBC Data Source Administrator

Click on the **System DSN** tab

Click Add

In the **Create New Data Source** Window

Select the Microsoft access driver (\*.mdb)

Click Finish

In the **ODBC Microsoft Access Setup** Window

Type “cdptcs-dsn” in the Data Source Name field

Leave Description field blank

Under Database, click the “Select” button

In the Directories column type C:\Program files\cdptcs

Select the **cdpdb.mdb** file

The Read Only and Exclusive boxes should be unchecked

Click “OK” for Select Database Window

Click “OK” in the ODBC Microsoft Access Setup Window

Click “OK” in the ODBC Data Source Administrator Window

Start> All Programs>CDPTCS

Change User Name to “Manager”

In the Console field, Select “01”

### **Windows XP Pro Notes:**

When installing on Windows XP machine Service Pack 2 and all Microsoft Updates must be installed on the machine.

For CDPTCS Setup Error Message “An error occurred while registering the file ‘C:\Windows\System32\msxml6.dll’. Select “Ignore”

### **Windows Vista Business Notes:**

Before running the CDPTCS installation program you must turn off the user account control

1. Use an Administrator Profile.
2. Turn off user account control

Go to Control panel > Click “User Accounts” > Click “User Accounts” again > Click “Turn User Accounts On or Off” > Uncheck “User Account Control (UAC) to help protect your computer” box and then restart computer

Start > Control Panel > (classic view) administrative tools.

All the Data Sources ODBC information is installed automatically.

### **Network Installation Notes:**

- I. Network Setup Installation
- II. The path to the server has to be shared
- III. Administrator must create a shared folder that can be seen on the network.
- IV. No drive mapping. And No IP addresses
- V. The CDPDB.mdb file must be put in the shared folder on the server.

Technical Support: (609) 599-2800

1/22/2008